

PC

Music

Phone

TV

Mobile

Remote

How to set up and demo ConnectLine™ with a Laptop



Checklist

What do I need?

- ✓ A pair of ConnectLine™ enabled hearing instruments. RITE styles are recommended so that domes can be replaced. Don't forget the batteries.
- ✓ Selection of domes.
- ✓ A Streamer (with power supply), and fully charged.
- ✓ Remember to pair the Streamer to the HIs in Genie, if you are not using a Demo Streamer.
- ✓ A complete ConnectLine™ TV kit (the variant sold in your market).
- ✓ A fully charged Bluetooth enabled mobile phone. If in doubt about version and compatibility – please visit the pairing guide on the Oticon website.
- ✓ Suitable video material with speech content. Either as a DVD (if your laptop has a DVD drive) or a DVD/Video file stored on the laptop.
- ✓ Laptop with audio output and video playback SW. Make sure that the Video player and sound settings provide audio out on the analogue output socket (a 3.5mm stereo mini jack – green).


Prepare Setup

1 Preparing the Demo

The preparation is pretty straight forward (as is the case with ConnectLine). Make sure that all electronics are charged, turned on, and that streamer, Mobile and TV adapter are pre-paired (see instructions below).

1 Pair Streamer and TV Adapter

Place the Streamer next to the TV adapter (less than 30 cm / 1 ft between devices).




2 Press the Streamer Connect button for 5 seconds until the blue light starts flashing.


3 Pairing is done when the blue light stops flashing rapidly. At the same time the STATUS indicator on the TV adapter and the Audio button on the Streamer will be orange.

1 Connecting the TV adapter

Connect the TV adapter to a wall socket. NB! using the USB cable to power the TV adapter directly from the Laptop might cause audible interference.



Connect the TV adapter to the green 3.5mm analogue audio output socket, using the cable and adapter included.



How to Demo

Pair Streamer and Mobile phone

- 1** Make sure that the Bluetooth function is turned on in the mobile phone. Find the menu or identify the activation required to start the pairing process. Refer to the manual of the mobile phone for detailed instruction on how to conduct the pairing.
- 2** Put Streamer into pairing mode by a very long press of 5 seconds on the Connect button until it flashes. Streamer is now in pairing mode indicated by the fast blinking blue light.
- 3** Activate the pairing feature on the mobile. The mobile is now searching for active Bluetooth devices in the area. When Streamer is found, the name 'Streamer 1.2' will appear in the display of the Mobile Phone. Go to the word 'Streamer 1.2' and press select or OK.
- 4** Streamer's code is: **0000** (four zeros). Upon entering the pin code, Streamer and the device are now paired.
- 5** After a successful pairing sequence, a symbol should appear in the mobile phone display (e.g. a headphone), indicating that a headset is now connected to the phone and that the phone can only be used via Streamer and hearing instruments.



How to Demo

- 1** Make sure that all practicalities are in place. PC powered on, mobile phone powered on with the Streamer connected (typically you will see a headset icon on your mobile).
- 2** NB! you will need a phone from which you can call your mobile phone. Make sure that you have access to one.
- 3** Make your customer wear the hearing instrument and the Streamer around his/her neck.
- 4** Explain that you will now demonstrate how the ConnectLine™ system enables seamless and intuitive switching from TV to Mobile phone and back. Don't forget to mention that the system works equally well with land line phones – thanks to the ConnectLine Phone adapter.
- 6** Start the video on the PC, and ask the dispenser to connect by pressing the audio button. Make sure that your customer can hear you and the "TV" signal.
- 5** Explain that you will now make a call. Call your own cell phone from any office phone available /outside the clinic/room.
- 7** Have a brief conversation and mention that all the user need to do to end the call and resume TV watching is to press the Phone button.
- 8** End the Demo – and discuss how you may help them give similar demos in their clinic.

People First

We believe that it takes more than technology and audiology to create the best hearing instruments. That's why we put the individual needs and wishes of people with hearing loss first in our development of new hearing care solutions.